



FDM ATR Polymers

- 500 polymers
- 100 plasticizers
- 170 monomers
- 160 trade name polymer additives
- 130 common and engineering plastics
- 80 rubber spectra named by functionality.
- 2 cm⁻¹ resolution
- Chemical and industrial names
- Many CAS numbers
- Many polymers have Mw, Mn and Tg values.
- Complements the FDM ATR Retail Adhesives & Sealants and the FDM Raman Polymers.


Poly (alpha-methylstyrene), 25014-31-7





Poly (isodecyl acrylate), Mw=60000, 5481-40-6



Poly (vinyl stearate), Mw=900000, 9003-95-6



Ethylene/vinyl acetate, Mw=100000, 24937-78-8



Fiveash Data Management, Inc.

608-236-9145

www.fdmspectra.com

Fax: 608-236-9170